
Indicadores Biológicos Biological Indicators

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
Parameters determined at time of manufacture according to EN ISO 11138-1: 2006, EN ISO 11138-2:2009 and IRAM 37102: 1999 (Parts 1 and 2). The showed values are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant

ATCC is registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

Industria Argentina - Made in Argentina

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Biological Indicators BT10 For Ethylene Oxide sterilization English

Composition

Each tube contains a population of spores of *Bacillus atrophaeus* ATCC 9372 soaked in a strip of paper. It also has a growth indicator media contained in the glass ampule.

Product description

Bionova® BT10 biological indicator is specifically designed for the monitoring of ethylene oxide (EO) sterilization processes. If the sterilization process was not successful, the indicator media will change to yellow after incubation at 37 ± 2 °C, thus indicating the presence of live *Bacillus atrophaeus* spores. If the sterilization process was successful, the indicator media will remain the original color after incubation. The final readout should be made after 48 hours of incubation at 37 ± 2 °C.

Warning!

Do not use BT10 biological indicator for monitoring steam sterilization cycles, dry heat, chemical vapor, radiation or other sterilization processes.
Do not re-use biological indicators.

Storage

Store it in a dark place at temperatures between 15-30 °C, 35-60 % relative humidity.
Do not freeze.
Do not store these biological indicator near sterilizing agents or other chemical products.

Directions for use

- Identify the Bionova® BT10 indicator tube by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.
- Pack the biological indicator along with materials to be sterilized in an appropriated package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for the sterilizing agent (EO). Generally a problematic area is the center of the load.
- Sterilize as usual.
- After the sterilization process has finished, you can:
 - Remove the biological indicator from the sterilization package for incubation prior to material aeration.
 - Aerate the package containing the biological indicator with the rest of the load and then remove the biological indicator for incubation.
- Check the chemical indicator on the label the biological indicator. A color change to green confirms that the biological indicator has been exposed to ethylene oxide. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.
- Crush the it ampule contained in the biological indicator and incubate at 37 ± 2 °C.
IMPORTANT: Use a non-sterilized biological indicator as positive control each time a processed indicator is incubated. The positive control ensure that correct incubation conditions were met.
- Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours at 37 ± 2 °C. Reading should be made at convenient intervals of 10 hours. A color change to yellow of the growth indicator media means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is made (the sterilization process was acceptable). The positive control indicator should show a yellow color change for the results to be valid.
Record the positive ones and discard them immediately as is shown below.
WARNING! Do not re-use the sterilizer until the biological indicator test results are negative (process indicator remain the original color).

Disposal

Dispose of used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.