Indicadores Biológicos Biological Indicators

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System Sensibility:

The sensibility of the system was determined as the difference between the positive indicators at 48 hs of incubation and the false negative (negative for the fluorescent readout but visually positive) in relation to those positive at 48hs of incubation.

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
Parameters determined at time of manufacture according to ISO 11138 (Parts 1 and 3) and IRAM 37102 (Parts 1 and 3) standards. The showed values are reproducible only under the same conditions under which it were determined.

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

Industria Argentina - Made in Argentina

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Biological Indicators BT110 For Ethylene Oxide sterilization

English

Composition

BT110 Bionova® Rapid Readout biological indicators for ethylene oxide sterilization processes have been designed for the fast and easy monitoring of ethylene oxide (EO) sterilization cycles. They consist of a plastic tube, a special filter, a spores carrier and a glass ampoule with culture medium. The tube has a plastic top with three holes and a barrier permeable to EO.

Rapid readout: 4 hours

The rapid readout must be carried out in the Bionova® Reader Incubator or in a similar reader incubator (read manufacturer's instructions for use). The reader-incubator must be capable of reading the fluorescence emission of the product, resulting from breaking a specific substrate on the spores carrier. Fluorescence is produced when the reader stimulates the carrier with UV light at 360nm. It is a direct process of the germination and growth of *Bacillus atrophaeus* spores which have survived the sterilization process (positive result).

Visual readout: 48 hours

If the sterilization process has not been successful, the purple culture medium will change to a greenish color and then to yellow, after incubation at 37±1 °C, thus showing the presence of live *Bacillus atrophaeus* ATCC 9372 spores. If the sterilization process is successful the culture medium will remain blue after the incubation process.

THE FINAL READOUT OF A NEGATIVE RESULT SHOULD BE MADE AFTER (4) FOUR HOURS OF INCUBATION. HOWEVER, IT IS ADVISABLE TO CONFIRM THE RESULT VISUALLY AFTER 48 HOURS OF INCUBATION.

Composition

Each tube contains a population of *Bacillus atrophaeus* ATCC 9372 spores soaked on a paper strip (spores carrier). It also has a culture media of blue color contained in the glass ampoule, over the filter and the carrier, on the base of the tube.

Warning!

Do not use BT110 biological indicators to control steam, dry heat, formaldehyde or other sterilization process. Do not reuse the biological indicators.

Storage

Store in dark place, at 15-30 °C temperature, 35-60 % relative humidity. Do not freeze.

Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

- Identify the BT110 Bionova® biological indicator by writing the sterilizer number (if you have more than one), load number, and processing date on the label.
- Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which a priori you consider most inaccessible for sterilizing agent (EO). Generally a problematic area is the center of the load and near the door.
- Sterilize as usual.
- After the sterilization process has finished, open the sterilizer door and remove the biological indicator from the package.
- Check that the chemical indicator printed in the label of the biological indicator has changed from purple to green. The color change indicates that the biological indicator has been exposed to the ethylene oxide. **IMPORTANT:** this color change does not evidence the process effectiveness to achieve sterility. If the chemical indicator color has not change, it is necessary to check the sterilization process.
- Crush the ampoule contained in the biological indicator and hold it by the top, beat the tube against a hard surface until the liquid wets the carrier. Incubate at 37±1 °C.
- IMPORTANT:** use a non-sterilized biological indicator as positive control every time a processed indicator is incubated. The positive control ensures that correct incubation conditions were met. Both the positive control indicator and the processed indicator should belong to the same batch. Before incubation, use the processed indicator to calibrate the reader, since it works through comparison. Do not calibrate the reader with an incubated indicator. The calibration should be carried out every time a biological indicator is incubated.
- Incubate the processed biological indicator and the indicator used as positive control for a maximum of 4 hours at 37±1°C for the rapid readout and during 48 hours for the color changes readout (visual result). Reading should be made at convenient intervals of 1 hour. The fluorescence detection done by the reader (stimulation 340-380 nm / emission 455-465 nm) or the growth indicator medium color change from blue to yellow means a failure in the sterilization process. If, after 4 hours (fluorescence detection) or after 48 hours (visual confirmation) no fluorescence or no color change is detected in the processed indicators, the result is negative (the sterilization process has been effective). The indicator used as a positive control should change from blue to yellow so as to consider the results valid.

Record the positive results and discard them immediately as it is explained below.

WARNING! Do not reuse the sterilizer until the biological indicator test result is negative (the processed indicator remains blue).

Disposal

Dispose used biological indicators according to your country's healthcare and safety regulations. Positive biological indicators can be autoclaved at 121° C for at least 15 minutes or at 132° C for 10 minutes in a gravity displacement sterilizer.