
Indicadores Biológicos Biological Indicators

Para la esterilización con VH_2O_2
For VH_2O_2 sterilization
Para a esterilização por VH_2O_2
Pour stérilisation avec VH_2O_2
Per sterilizzazione con VH_2O_2

BT90

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Certificado de calidad Quality Certification Bionova® BT90

Esterilización por VH_2O_2 / VH_2O_2 sterilization
Geobacillus stearothermophilus ATCC 7953

Población / Population _____ UFC / CFU

Valor D / D-value _____ min.

D-value (2 mg/L H_2O_2 Vapor, 45 °C)

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
Parameters determined at time of manufacture according to EN ISO 11138-1: 2006 and IRAM 37102-1: 1999. The showed values are reproducible only under the same conditions under which it were determined.

ISO and USP Compliant
ATCC is registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

Industria Argentina - Made in Argentina

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English

Description

Bionova® F1-BT90 Rapid Readout Colorimetric System Biological Indicators for Plasma have been designed for the easy evaluation of sterilization processes. The system contains a tube with a conic base and a special filter that keeps the glass at the top of the tube. An easy readout can be carried out using the reference color table.

If the sterilization process was not successful, the indicator media will change from purple to grey-green and finally to yellow after incubation at 60 ± 2 °C, thus indicating the presence of live *Geobacillus stearothermophilus* ATCC 7953 spores.

If the sterilization process was successful, the indicator media will remain purple after incubation. The final readout should be made after 8 hours of incubation at 60 ± 2 °C.

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a carrier. It also has a purple growth indicator media contained in the glass ampule on a plastic filter.

Warning!

Do not use F1-BT90 biological indicators for EO, Steam, Dry Heat, Formaldehyde or other sterilization process. Do not re-use the biological indicators.

Storage

Store it in a dark place and at temperatures between 15 - 30 °C, 35 - 60 % relative humidity.

Do not freeze.

Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use

1. Identify the Bionova® F1-BT90 indicator by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.

2. Pack the biological indicator along with the materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for the sterilizing agent. Generally a problematic area is the center of the load and near the door.

3. Sterilize as usual.

4. After the sterilization process has finished, open the sterilizer door and remove the biological indicator from the package. **WARNING!** Do not crush or handle the glass ampule too roughly.

5. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to VH_2O_2 plasma agent. **IMPORTANT:** this color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

6. Use the clamp to press the tube at the level of the tag and crush the ampoule hidden above the plastic filter and incubate the biological indicator at 60 ± 2 °C.

IMPORTANT: Use a non-sterilized biological indicator as a positive control each time a processed biological indicator is incubated. The positive control ensures that correct incubation conditions were met.

7. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 8 hours at 60 ± 2 °C. Reading should be made at convenient intervals of 1 hour.

A color change of the growth indicator media from purple to grey-green means a sterilization process failure has occurred. If after 8 hours there is no color change in the processed indicators, a final negative result is made (the sterilization process was acceptable). The positive control indicator should show a purple to grey-green color change for the result to be valid. The final color should be yellow.

DETECTION METHOD: THE SYSTEM PERMITS AN EARLY DETECTION OF THE POSITIVES BY THE GRADUAL COLOR CHANGE OF THE CULTURE MEDIUM.

The readout must be carried out using the reference color table, under the white light source. The yellow lights may hinder the readout. The table permits the easy detection of positive biological indicators. The use of the reference color table allows a final lecture within the first 8 hours of incubation with a 97 % of sensibility.

Record the results and discard immediately as it is shown below. **WARNING!** Do not use the sterilizer until the biological indicator test results are negative.

Disposal

Dispose of used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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