Indicadores Biológicos Biological Indicators

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Spores ampoules
For Steam sterilization

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
Parameters determined at time of manufacture according to ISO 11138 (Parts 1 and 3) and IRAM 37102 (Parts 1 and 3) standards. The showed values are reproducible only under the same conditions under which it were determined.

ISO and USP Compliant.

ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Quality Assurance Director

Industria Argentina - Made in Argentina

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Biological Indicators For Steam sterilization

English

Composition

Each Bionova® BT22 self-contained biological indicator contains a *Geobacillus stearothermophilus* ATCC 7953 spores population in 0.3 ml of culture media encased in a glass ampule (borosilicate type 1) which has a diameter of 7.0 mm, wall thickness of 0.2 - 0.3 mm and a length of 25.0 mm. The culture media will change to yellow to show signs of growth. Also, every ten units, there is one negative control belonging to the same batch as the culture medium but without spores.

Product description

Bionova® BT22 biological indicator is specifically designed for the monitoring of steam sterilization processes at 121 - 134 °C.

If sterilization process was not successful, the indicator media will change to yellow after incubation at 56±2 °C, thus indicating the presence of living *Geobacillus stearothermophilus* spores.

If the sterilization process was successful, the indicator media will remain the original color after incubation. The final readout should be made after 48 hours of incubation at 56±2 °C.

The negative control can be used in two ways:

1. It can go with the spore ampoules during the process. In this way, it will verify that any color change taking place after the sterilization process is a result of it and not of the germination and growth of spores.

2. It can be incubated without being exposed to the sterilization process, together with the positive control (a spore ampoule not exposed to the sterilization process) and the sterilized ampoules. In this way, it is used as a color pattern to be compared with the positive and negative results.

Warning!

Do not use BT22 biological indicator for monitoring ethylene oxide, dry heat, chemical steam, radiation or other sterilization processes. Do not re-use the biological indicators.

Storage

Store it in a dark place and at temperatures between 4 - 8 °C.

Do not freeze.

Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use

1. Take the Bionova® BT22 ampoules out of the refrigerator and leave them at room temperature for at least one hour.

2. Place the Bionova® BT22 ampoules inside a container with the liquid being sterilized or inside a similar container with liquid. This container should be of the same size and have the same volume as the one containing the liquid being sterilized.

IMPORTANT: You should use at least ten indicators to monitor each sterilization cycle, together with one negative control.

3. Sterilize as usual.

4. Once the sterilization process has finished, open the sterilizer door, wait five minutes and then remove the biological indicator from the container. **WARNING!** Wear safety glasses and gloves when removing the BT22 biological indicator from the sterilized container.

5. Incubate the ampoules at 56±2 °C.

IMPORTANT: Use a non-sterilized biological indicator as positive control each time a processed biological indicator is incubated. The positive control ensures that correct incubation conditions were met. 6. Incubate the processed biological indicator, the negative control and the indicator used as positive control for a maximum of 48 hours at 56±2 °C. Reading should be made at convenient intervals of 10 hours.

A color change to yellow of the growth indicator media means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is made (the sterilization process was acceptable). The positive control indicator should change to yellow for the result to be valid. Record the positive ones and discard them immediately as it is shown below.

WARNING! Do not use the sterilizer until the biological indicator test results are negative (process indicator reads the original color).

Disposal

Dispose of used Bionova® BT22 biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.