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Class 5 Indicator according to EN ISO 11140-1:2009

Product Description

IT26 indicator integrator is three phased sequential indicator, developed to react in steam sterilization processes. The pink ink was fixed to change to black when the population of spores is completely dead, indicating **SAFE** condition is reached.

Level 1 – **Blue: Exposure level**, It turns to black to indicate the exposure to the process but it doesn't assure that the conditions indicated in the emulador have been met (**UNSAFE** condition).

Level 2 – **Pink: Integration level**, it turns to black thus indicating that conditions indicated in the emulador have been met. These conditions are calibrated with the kill time of a *Geobacillus stearothermophilus* 10⁶ spores population, calculated in BIER (Biological Indicator Evaluator Resistometer). This is called **SAFE** condition.

Level 3 – **Green: Over-integration level**, turns to black when integration conditions indicated in the product have been exceeded. This is called **OPTIMUM** condition.

Characteristics

Consists of a paper strip of 60 mm x 25 mm with brilliant OPP lamination on one side. Free of lead and toxic heavy metals.

Indications for Use

Use the indicator integrator for monitoring Steam sterilization processes. **SAFE** condition is reached in steam sterilization processes as described in the quality certificate.

Precautions

Do not expose this product to EO, dry heat, or any chemical or physical agent sterilization processes other than Steam. Do not reuse the sterilizer until the integrator's color reaches or passes the **SAFE** condition.

Instructions of Use

- Place an INTEGRON® IT26 indicator integrator within each package, bag or tray that will be sterilized by steam.
- In rigid containers place an INTEGRON® IT26 indicator integrator in every corner, or at least two diagonally opposite corners, to ensure that the steam has reached the most inaccessible areas.
- After completing the process, analyze the results. If the **SAFE** integration condition has been reached or exceeded, the sterilization process was successful.

Storage and Shelf life

- Indicators should be stored in dry places, protect from light at temperatures between 15 - 30 °C, 35 - 60% relative humidity. Do not wet. Do not store close to sterilizing agent.
- Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. The expiration date as the batch of the product are specified on the package label and in each unit of product. Do not use indicators after their expiration date.
- After use, the indicator will not change in 1 year if stored in temperature and relative humidity previously indicated.

Industria Argentina - Made in Argentina

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