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CHEMDYE®
BDA5
Español

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Bowie - Dick Test Sheet
BDA5
English

Product description

The Chemdye® Bowie-Dick test sheets are A5 sized printed with a chemical indicator in a concentric pattern. The chemical indicator turns from yellow color to dark brown / black color indicating adequate air removal. An air removal fault is displayed as a lighter colored area near the center of the otherwise dark brown/black test sheet.

Indication for use

The Chemdye® test sheets are designed to be used as the indicator placed inside facility-assembled standard Bowie-Dick test packs for testing air removal efficiency of 270 - 273 °F (132 - 134 °C) vacuum assisted steam sterilizers. The test pack should be placed in the horizontal position on the lowest shelf in an already warmed, empty sterilizer. The sterilizer should run for 3.5 - 4.0 minutes at 270 - 273 °F (132 - 134 °C). Sterilizers which show non-uniform color development should be shut down and scheduled for service.

Warning!

Do not use the Chemdye® test sheets for air removal monitoring of gravity displacement steam, dry heat, ethylene oxide, or other low temperature sterilization processes.

Storage

Store in a dark place and at temperatures between 15 - 30 °C, 30 - 60 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. Perform the air removal test daily and after sterilizer repairs.
2. Prepare a test pack of absorbent surgical towels. The towels should be freshly laundered, but not ironed. The towels must be folded to a size not smaller than 6.7 inches (17 cm) in one direction and 8.7 inches (22 cm) in the other direction, and placed one above another. The height of the test pack must be 10 to 11 inches (25 to 28 cm).
3. Place a single Chemdye® Bowie-Dick test sheet in the center of the test pack.
4. Place the towels and test sheet in a single, loosely applied wrap and tape securely.
5. Place pack in a horizontal position in the front/bottom of an empty vacuum assisted steam sterilizer.
6. Run the sterilizer for 3.5 - 4.0 minutes at 270 - 273 °F (132 - 134 °C).
7. After completion of the cycle, retrieve the Chemdye® Bowie-Dick test sheet and examine:

.Satisfactory Test: This is indicated by a test sheet that shows uniform dark brown / black color development. Such a test indicates rapid steam penetration and adequate air removal.

.Unsatisfactory Test: This is indicated by a light colored area in the middle of the dark brown / black test sheet. This light colored area is caused by an air pocket and indicates inadequate or incomplete air removal.

A failure with a Bowie-Dick test should be reported immediately to the supervisor. The sterilizer should not be used until examined for possible malfunction.

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