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English

Description

Chemdye® Bowie-Dick Test Packs detect air leaks, inadequate air removal and steam penetration for vacuum-assisted steam sterilizers. Bowie-Dick Test Packs should be used daily, after a sterilizer is installed, or relocated, after a sterilizer malfunction, after sterilization process failures, and after any major repairs of the sterilizer. The Chemdye® Bowie-Dick Test Pack consists of one chemical indicator sheet (12,5 x 12,5 cm) positioned inside porous materials and sealed inside a disposable outer wrap, which changes from yellow to black when is processed.

Any heterogeneous color change, such as the center of the test sheet's color being paler than the edges or a non-uniform color change indicates that there was an air pocket present during the cycle showing a malfunction in the process.

The usage of the Chemdye® BD TP must be carried on every day the sterilizer is used before the first process load. A shortened cycle (i.e. a cycle omitting the postvacuum drying phase) should be run first to properly heat up the sterilizer.

Instructions for use

1. Place one Bowie-Dick Test Pack in the horizontal position with the label facing up in an empty pre-vacuum sterilizer chamber, in the bottom/front of the sterilizer rack, and over the drain.
2. Run the sterilizer for not more than 17,0 minutes at 121 °C. If the test pack is exposed for more than 17,0 minutes, it will produce invalid results.

Interpretation of results

After completing the cycle, remove and examine the chemical indicator sheet.

Satisfactory Test: It is indicated by a test sheet that shows uniform black color in the extension of the indicator. This result points a rapid steam penetration and adequate air removal of the sterilizer chamber.

Unsatisfactory Test: It is caused by an air pocket within the chamber and indicates incomplete air removal from the chamber, poor quality wet steam, superheated steam and/or non-condensable gas in steam supply.

If a non-uniform color change is noted, it should be reported immediately to the supervisor and that sterilizer should not be used until it is examined. Put in quarantine the previous loads sterilized in that machine until the it's certain that the machine has not failed before.


Warning!

Do not use a Bowie-Dick Test Pack in a gravity displacement steam sterilization cycles (where there is no prevacuum phase to check) or in other kinds of sterilization processes.

Storage

Store Chemdye® Bowie-Dick Test Packs in a dark place, 30 - 60 % RH at 15 - 30 °C.

The Bowie-Dick Test Pack is manufactured in accordance to the EN ISO 11140-1: 2009 and EN ISO 11140-4: 2007.

Industria Argentina - Made in Argentina

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