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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Parameters determined at time of manufacture according to ISO 11138-1: 2006, ISO 11138-3: 2006, ISO 11138-5: 2006 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Composition

Each Bionova® BT50 envelope consists of a paper strip soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product description

The envelope containing Bionova® BT50 Spore Strip is specifically designed to control Steam 121 °C-134 °C (Steam) and Formaldehyde (Form) sterilization processes. Spore strips should be incubated in Bionova® MC20 or MC1020 Growth Medium tube or in appropriate Growth medium for *Geobacillus stearothermophilus* ATCC 7953 growth. If the sterilization process was not successful, MC20 or MC1020 Growth Medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores on the strips.

If the sterilization process was successful, MC20 or MC1020 Indicating Medium will remain its original color. A final readout should therefore be performed after 24 hours of growth medium incubation at 60±2 °C for Steam sterilization control processes, or 48 hours at 60±2 °C for Formaldehyde sterilization control processes.

Warning!

Do not use Bionova® BT50 Spore Strip envelopes to control EO, Dry Heat sterilization, Radiation or other sterilization processes different from those for which were designed. Do not reuse spore strip envelopes.

Storage

Store in a dark place and at temperatures between 15-35 °C, 35-60 % Relative Humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. Identify Bionova® BT50 envelope by writing on its back the sterilizer number (in case of having more than one), load number and processing date.

2. Pack it in an appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the envelope in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.

3. Sterilize as usual.

4. After sterilization process has finished, remove Bionova® BT50 envelope from the package for processing and incubation.

5. Tear up the envelope at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer the spore strip with a sterile clamp to MC20 or MC1020 Growth Medium tube or any other appropriate growth medium. **IMPORTANT!** Use latex gloves and chin strap when transferring spore strips from the envelope to growth medium tube. **WARNING!** Avoid spore strip contact with any kind of surface including the outer part of Growth Medium tube and of Bionova® BT50 envelope.

IMPORTANT: Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

6. Incubate processed spore strips along with the strip used as a positive control at 60±2 °C during a maximum of 24 hours for Steam or 48 hours for Formaldehyde. It is advisable to make observations every 10 hours.

Color change of MC20 or MC1020 growth medium to yellow means a failure of the sterilization process. If after 24 hours (Steam) or 48 hours (Form) no color change is visible on MC20 or MC1020 growth medium, the result is negative (i.e. the sterilization process was successful). MC20 or MC1020 Indicating Medium used for incubating the spores strips which have not been sterilized (positive control) should turn yellow in order for the results to be valid. Record the positive ones and dismiss them immediately as indicated below. **WARNING!** Do not reuse sterilizer until spore strip growth result is negative (MC20 or MC1020 Growth Medium containing processed spore strip remains its original color).

Disposal

Discard Growth Medium tubes after use to analyse Bionova® BT50 Spore Strips growth according to your country's healthcare and safety regulations. Positive MC20 or MC1020 Growth Medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.