MC1030-2

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MC1030-2

Culture medium with pH indicator.

For growth of *Bacillus atrophaeus* ATCC 9372 and *Bacillus subtilis* DSM 5230 bacterial spores.

Usage

To test growth of *Bacillus atrophaeus* ATCC 9372 bacterial spores used for monitoring Dry Heat and EO sterilization processes and *Bacillus subtilis* DSM 5230 used for monitoring Steam at low temperature sterilization processes.

Applicable Regulation

ISO11138-1:2006, ISO11138-2:2006, ISO11138-3:2006 and ISO11138-4:2006; IRAM 37102-1:1999 and IRAM 37102-2:1999.

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

Characteristics

Sterile growth medium is specially formulated for use with our BIONOVA® products (also suitable for other brand products).

Growth medium color: blue

Plastic tube: 68 mm high, 10 mm diameter.

Medium volume: 2 ml

For germination and growth of *Bacillus atrophaeus* ATCC 9372 and *Bacillus subtilis* DSM 5230 bacterial spores.

Incubation conditions:

48 hours (maximum). 37± 2 °C.

Medium turns to yellow when spores grow.

Environmental conditions during manufacture

15-30 °C room temperatures, RH 30-80 %. Sterility conditions are necessary during plastic tube filling with sterile medium, performed in laminar flow.

Storage conditions

T = 4-20 °C, RH 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed.

Transport in closed and reinforced boxes in order to avoid damages.

Transportation of this product does not represent a risk for health.

Shelf-life

24 months.

Packing

50 units per box.

Information only on product's label.

Labelling

On product: Label 17 mm x 33 mm printed in black. Product code, batch number and bacterial strain.

On product's box: product code and description, presentation, operating instructions, batch number, bacterial strain and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

It is advisable to incubate at 37 °C in BIONOVA® IC10/20 incubator.

Note: When necessary, limits of each product's parameters should be specified.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not freeze.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.