# BT24S/6

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# BT24S/6

## Spore suspension

For direct inoculation of products or preparation of customized biological indicators for monitoring Steam sterilization of liquids at low temperatures (110 °C – 121 °C).

## Usage

For monitoring low temperature (110 – 121 °C) Steam sterilization of liquids.

## Applicable Regulation

ISO 11138-1:2006 and ISO 11138-3:2006; IRAM 37102-1:1999 and IRAM 37102-3:1999.

## Classification

Class 1, according to risk.

## Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Technology) PM 1614-1.

## Characteristics

Glass bottle: 18 mm diameter, 0.3 mm thick wall, 60.7 mm high.

Plastic cap: 16 mm external diameter, 9.4 mm high.

Population:  $10^5$  *Bacillus subtilis* DSM 5230 ATCC 35021 spores per 0.1 ml ( $10^6$  spores per ml), prepared in 40% ethanol.

D-Value<sub>STEAM</sub> is provided. Conditions: 115 °C. Another D Value is calculated at another temperature for Z value calculation.

Z Value<sub>STEAM</sub> is provided.

Upon completion of the sterilization cycle, the inoculated product should be tested in a suitable culture medium at 37 °C. We suggest the use of culture medium Bionova® MC1020 or MC1030 to final reading in 48 hours. For conventional culture media like TSB, incubation for 7 days is recommended.

## Environmental conditions during manufacture

Refrigeration during inoculation. Sterility conditions are necessary only during the filling and sealing process, carried out under laminar flow conditions.

## Storage conditions

For short-term use, store at 2–8 °C in refrigerator. Freeze to store for long periods (-10 to -25 °C). RH: 35-60%. Keep out of direct light.

## Transportation conditions

Storage conditions should be strictly followed.

Transport in closed reinforced boxes to avoid damage. This product does not represent a risk for health. Bottles can be stored at temperature below 30 °C for up to 15 days during transportation.

## Shelf-life

24 months.

**Packing**

1 bottle with 10 ml of spore suspension within a plastic bag.

**Labelling**

On the bottle: Product code and description, presentation and spore concentration, directions for use, batch number, manufacture and expiration date and manufacturer information.

**Possible target markets**

Healthcare and Industry

**Other important information (optional)**

It is recommended to incubate at 37 °C in BIONOVA® IC10/20 incubator.

Read product's directions for use thoroughly before use.

**Note:** When necessary, limits of each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be sterilized in gravity air displacement steam sterilizers at 121°C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

Note: Spores can be incubated in culture media BIONOVA® MC1030 or MC1020 culture media or other culture medium made up of soybean casein and yeast extract.